
An In Vitro Study of the Treatment of Implant Surfaces With Different Instruments

Reiner Mengel, Dr med dent, DDS*/Christian-Eiben Buns, Dr med dent, DDS*/

Claudia Mengel, Dr med dent, DDS*/

Lavin Flores-de-Jacoby, Prof Dr med dent, PhD, DDS*

The aim of this in vitro study using implants and abutments was to examine traces left by various cleaning instruments and to determine the quantity of substance removal. The implants and relevant abutments were Screw-Vent implants (Dentsply), titanium plasma-coated full-screw implants (Straumann), and standard Brånemark implants (Nobel Biocare). The instruments used were titanium curettes, Gracey curettes, plastic curettes, rubber cups with Zircate prophylaxis paste, the Cavitron Jet ultrasonic scaler with universal inserts and airpolishing nozzles with Prophyl-Jet cleaning powder, and the Denson sonic scaler with SofTip disposable prophyl tips and universal tips. Two test fields (2 × 2 mm) on the respective implant head and abutment underwent standard treatment with each instrument. With the full-screw implants, an additional two test fields were examined from the transition of the titanium plasma coating to the implant head. Untreated implants and abutment surfaces served as controls. The treated and untreated surfaces were examined by scanning electron microscopy for work traces and by optical laser profilometry for substance removal. The results showed that the Cavitron Jet ultrasonic unit, the Denson sonic scaler with universal tips, and the Gracey curette left pronounced traces and induced high substance removal. The titanium curette and the Denson sonic scaler with SofTip disposable prophyl tips left virtually no traces and removed very little substance. The rubber cup, the plastic curette, and the Cavitron Jet airpolishing system left the implant surfaces unchanged. All instruments apart from the rubber cup and Cavitron Jet airpolishing system left pronounced traces at the transition of the implant head to the titanium plasma coating of the full-screw implants. The Cavitron Jet airpolishing system, the rubber cup, the plastic curette, and, with some reservations, the Denson sonic scaler with SofTip plastic fittings and the titanium curette, appear to be suitable for cleaning implant surfaces.

(INT J ORAL MAXILLOFAC IMPLANTS 1998;13:91-96)

Key words: implant surface, profilometry, scanning electron microscopy

The main etiologic factor in periodontal disease is plaque, as the bacteria present in it are responsible for inflammation of the periodontal tissue.¹ It is therefore crucial to remove plaque at an early stage. Although this task is primarily the responsibility of the patient, regular plaque controls should be carried out by the dentist, with any plaque or calculus

residues being removed. The methods that have been proven successful in such applications are manual,² mechanically rotating² or oscillating³ instruments, ultrasonic equipment, airpolishing nozzles with cleaning powder, and sonic scalers.

The main etiologic factor in peri-implant inflammation is also plaque or its bacterial components.⁴ The peri-implant area seems to be even more susceptible than the periodontium⁵ to bacteria, indicating that early plaque removal is even more essential in implant-treated patients.⁶ Although the patient can remove plaque with standard prophylactic agents (eg, toothbrush, dental floss, interdental brush), no instruments are available without referral to the dentist for professional implant cleaning.

*Department of Periodontology, Marburg University, Marburg, Germany.

Reprint requests: Dr Reiner Mengel, Dental School of Medicine, Department of Periodontology, Georg-Voigt-Str. 3, D-35033 Marburg, Germany.




Fig 1 Device used to carry out the reproducibility of the experiments.

The aim of this *in vitro* study using implants and abutments was to examine by scanning electron microscopy the work traces left by various instruments, and to determine by optical laser profilometry the quantity of substance removal.

Materials and Methods

The following implants and abutments were tested:

1. Screw-Vent implant (two-part titanium screw implant) with Screw-Vent abutment (Dentsply, York, PA).
2. ITI full-screw implant (one-part titanium plasma-coated titanium screw implant) (Straumann, Waldenburg, Switzerland)
3. Standard Brånemark implant (two-part titanium screw implant) with standard abutment (Nobel Biocare, Göteborg, Sweden)

The implants and abutments used were new standard products, which were deployed without the addition of a plaque-like substance. This study focused on pure removal of the titanium surfaces to permit comparison of the material loss after treatment with various instruments.

The following instruments were used:

1. Titanium curette (Deppeler SA, Rolle, Switzerland)
2. Gracey curette (Hu-Friedy, Chicago, IL)
3. Plastic curette (Nobel Biocare)
4. Rubber cup with Zircate prophylaxis paste (Dentsply)
5. Cavitron Jet ultrasonic scaler with universal insert (Dentsply)
6. Cavitron Jet airpolishing nozzle with Prophy-Jet cleaning powder (Dentsply)
7. Densononic sonic scaler with SofTip disposable prophylaxis tip (Dentsply)

8. Densononic sonic scaler with universal tip (Dentsply)

The instruments were used following the recommendations of the respective manufacturers. With each of the three different curettes and the rubber cup, five procedures were performed for 20 seconds at a speed of 5000 rpm. The universal insert of the Cavitron Jet ultrasonic unit was used at 30,000 Hz, and the Densononic sonic scaler with the universal and disposable tip at 3 bars pressure and an operating frequency of 6300 Hz for 20 seconds, respectively. The Cavitron Jet airpolishing nozzle with Prophy-Jet cleaning powder was operated from a distance of 2 to 3 mm for 20 seconds at a pressure of 3 bars. For handling purposes, the curettes, the rubber cup, the Densononic sonic scaler, and the Cavitron ultrasonic unit were firmly clamped in an arm balanced with weights to 0.2 N. The implants and abutments were fixed with a special-purpose bonding agent on a horizontally movable slide to ensure standardized treatment.

Two test fields (each 2×2 mm) on each implant head and at the center of each abutment were treated with each instrument (Fig 1). In the case of the ITI full-screw implant, two test fields at the transition of the titanium plasma coating to the implant head were treated in addition to the implant head.

Untreated implant and abutment surfaces in the immediate vicinity of the test fields served as control fields. A total of 24 implants and abutments were evaluated.

Following treatment, one test and one control field per implant or abutment were sputter-coated with gold in an S150 sputter coater (Edwards, Frankfurt, Germany), and then examined for work traces under an ISI-SX-30 scanning electron microscope (SEM) (LEO-Elektronenmikroskopie, Oberkochen, Germany). The quality of the work traces on the treated surfaces was assessed in four grades (none, slight, moderate, and pronounced) by two independent testers, with the untreated surfaces serving as controls. Prior to SEM, substance removal from the other test field was determined four times with a mechanical profilometer (Hommeltester T6D, Hommelwerke, Schweningen, Germany). This measuring and registration unit, which records surface form deviations of the first to fourth order as per German standard (DIN) 4760, operates with an inductive sensor with a fine diamond tip ($2\text{-}\mu\text{m}$ radius). The surface probe was drawn over the surface to be marked by means of a precision feed unit in selectable scanning paths and at constant scanning speeds, and the scanned surface profile was converted directly into electric voltage proportionate in size and polarity. The carrier frequency principle

Table 1 Average Work Traces and Substance Removal on Implants and Abutments After Treatment With Different Instruments

Instrument	Work traces*	Substance removal (μm)		
		Brånemark	Screw-Vent	ITI
Cavitron with airpolishing nozzle	None (0)	0.05 ± 0.01	0.04 ± 0.01	0.07 ± 0.02
Rubber cup	None (0)	0.1 ± 0.02	0.1 ± 0.03	0.1 ± 0.03
Plastic curette	None (0)	0.1 ± 0.02	0.1 ± 0.02	0.1 ± 0.01
Densonc with SofTip	Slight (1)	0.15 ± 0.03	0.2 ± 0.04	0.2 ± 0.04
Titanium curette	Slight (1)	0.2 ± 0.05	0.2 ± 0.04	0.2 ± 0.04
Gracey curette	Pronounced (3)	0.5 ± 0.06	0.6 ± 0.06	0.5 ± 0.05
Densonc with universal tip	Pronounced (3)	0.8 ± 0.07	0.9 ± 0.10	0.8 ± 0.09
Cavitron with universal tip	Pronounced (3)	1.1 ± 0.10	1.2 ± 0.10	1.2 ± 0.20

*Roughness and loss of implant substance were evaluated by a modified Roughness Loss of Tooth Substance Index (Lie and Leknes18) according to the following criteria: 0 = smooth and even implant surface without visible marks from instrumentation and with no loss of implant substance; 1 = slightly roughened or corrugated implant areas; 2 = moderately roughened or corrugated implant areas; 3 = pronounced roughened or corrugated implant areas with a considerable number of lesions from instrumentation.




Fig 2 Abutment (Brånemark) treated with a Gracey curette. Pronounced work traces in handling direction and high substance removal (original magnification $\times 201$).




Fig 3 Abutment (Brånemark) treated with a Cavitron Jet ultrasonic scaler with universal insert. Pronounced single work trace on an otherwise smooth surface and high substance removal (original magnification $\times 200$).

applied here permits both static and rapidly responding dynamic measurements. The high-precision measuring amplifier is capable of registering and evaluating minute readings of up to 10^{-5} mm. The analog measurement signal is amplified up to 10^{-5} fold in the measuring amplifier, guided over a wave filter, evaluated in the computing circuit, and held on call in six measured-value memories. In the present study, the profile height (PH) yielded by the distance between the upper and lower reference line in relation to the scanned path was calculated according to DIN 4771.

Results

The Gracey curette, the Cavitron Jet ultrasonic scaler, and the Densonc sonic scaler with universal tip left moderate to pronounced work traces and caused an average substance removal of $0.83 \mu\text{m}$

(Table 1, Figs 2 to 4). The average substance removal measured after treatment refers to the loss of substance on the titanium surface itself. The titanium curette and the Densonc sonic scaler with SofTip disposable prophylaxis tip left only slight working traces and caused an average substance removal of $0.19 \mu\text{m}$ (Table 1, Figs 5 and 6). The rubber cup, the plastic curette, and the Cavitron Jet airpolishing system caused no visible change to the implant surfaces, recording an average substance removal of $0.09 \mu\text{m}$ (Table 1, Figs 7 to 9). All instruments, apart from the rubber cup and the Cavitron Jet airpolishing system, left pronounced traces at the transition of the implant head to the titanium plasma coating of the full-screw implants (Figs 10 and 11).

Following treatment with all instruments, no difference was recorded between the different implant systems and between the implants and abutments with respect to work traces and substance removal.




Fig 4 Abutment (Screw-Vent) treated with a Denasonic sonic scaler with a universal tip. Inhomogeneous rough surface with pronounced work traces and very high substance removal (original magnification $\times 200$).




Fig 5 Abutment (Screw-Vent) treated with a titanium curette. Homogeneous surface with slight work traces in handling direction and very low substance removal (original magnification $\times 200$).




Fig 6 Abutment (Screw-Vent) treated with a Denasonic sonic scaler with SoftTip disposable prophylaxis tip. Uniform surface with slight work traces and very low substance removal. The deposits on the surface are pieces of the disposable tip (original magnification $\times 200$).




Fig 7 Abutment (Bränemark) treated with a rubber cup and Zircate prophylaxis paste. No visible surface changes (original magnification $\times 200$).




Fig 8 Implant head (ITI) treated with a plastic curette. No visible surface changes (original magnification $\times 200$).




Fig 9 Implant head (ITI) treated with a Cavitron Jet airpolishing nozzle with Prophylaxis-Jet cleaning power. No visible surface changes (original magnification $\times 200$).




Fig 10 Untreated transition from titanium plasma coating to implant head (ITI) (original magnification $\times 200$).




Fig 11 Treated transition from titanium plasma coating to implant head (ITI). Pronounced work traces resulting from treatment with a Gracey curette (original magnification $\times 200$).

Discussion

Because of their high free-surface energy, titanium implants appear to accumulate more plaque than natural teeth,⁷⁻⁹ with up to 25 times more bacteria adhering to rough than to smooth implant surfaces.⁷ Any damage to the surfaces of the implant or abutment induces a change in the chemical oxide layer, which may ultimately result in increased corrosion. This process impairs the adhesion of fibroblasts, jeopardizing the biocompatibility of the implant.¹⁰ These findings have led to a demand for plaque and calculus removal from implants only with instruments that do not cause surface damage. In particular, Gracey curettes,¹⁰⁻¹⁵ sonic scalers with universal tips,^{12,15,16} and ultrasonic scalers with universal inserts^{15,16} should be avoided, since they leave pronounced work traces on the implant surface. Although better results were achieved with Teflon-coated inserts for sonic and ultrasonic scalers, they too should be used with caution.¹⁵ This warning also applies to the use of titanium curettes, which may also leave slight work traces on the implant surface.^{10,16} In contrast, virtually no work traces were left by plastic curettes,¹⁰⁻¹⁵ airpolishing nozzles with cleaning powder,^{12,13,16,17} and rubber cups.^{12-14,16}

The present results are comparable to those of other studies. In addition, the high figures recorded for substance removal by the metallic instruments demonstrate that they are not suitable for cleaning implants. Although titanium curettes and the Densonic sonic scaler with SofTip disposable prophylaxis tips are more suitable, and despite the fact that they leave only slight work traces, substance removal at the implant surfaces is more substantial. As this study involved using each of the instruments only once, the amount of substance lost by using titanium curettes

and the Densonic sonic scaler with SofTip disposable prophylaxis tips within the scope of a multiyear recall can only be estimated. These instruments should be used only for removing coarse deposits, but not for removing residual calculus and plaque, which should be done with a plastic curette and/or a rubber cup, as these instruments cause no visible changes to the implant surface. Although no work traces or substance removal were caused by the Cavitron Jet ultrasonic scaler with Prophy-Jet cleaning powder, it cannot be recommended without reservation, since the risk of trauma to the peri-implant tissue cannot be excluded.

Summary

The rubber cup, plastic curette, and Cavitron Jet airpolishing system can be used for supragingival removal of calculus and plaque on implant surfaces without the risk of damage. Further investigations are needed to determine which instruments are suitable for cleaning the subgingival peri-implant area.

References

1. Listgarten MA. Pathogenesis of periodontitis. *J Clin Periodontol* 1986;13:418-425.
2. Mengel R, Stelzel M, Flores-de-Jacoby L. An in vitro study of various instruments for root planing. *Int J Periodont Rest Dent* 1997;17:593-599.
3. Mengel R, Buns C, Stelzel M, Flores-de-Jacoby L. An in vitro study of oscillating instruments for root planing. *J Clin Periodontol* 1994;27:573-578.
4. Berglundh T, Lindhe J, Marinello C, Ericsson I, Liljenberg B. Soft tissue reactions to de novo plaque formation on implants and teeth. An experimental study in the dog. *Clin Oral Implants Res* 1992;3:1-8.
5. Ericsson I, Berglundh T, Marinello CP, Liljenberg B, Lindhe J. Long standing plaque and gingivitis at implants and teeth in the dog. *Clin Oral Implants Res* 1992;3:99-103.

6. Quirynen M, van Steenberghe D. Bacterial adhesion to oral implants and assessment of attachment and marginal bone level. *Dtsch Zahnärztl Z* 1993;48:158-160.
7. Quirynen M, Bollen CML. The influence of surface roughness and surface-free energy on supra- and subgingival plaque formation in man. A review of the literature. *J Clin Periodontol* 1995;22:1-14.
8. Quirynen M, Van Der Mei C, Bollen CML, Schotte A, Marechal M, Doornbusch GI, et al. An in vivo study of the influence of the surface roughness of implants on the microbiology of supra- and subgingival plaque. *J Dent Res* 1993;72:1304-1309.
9. Quirynen M, Van Der Mei C, Bollen CML, Van Den Bossche LH, Doornbusch GI, van Steenberghe D, et al. The influence of surface-free energy on supragingival plaque microbiology. An in vivo study on implants. *J Periodontol* 1994;65:162-167.
10. Dmytryk JJ, Fox SC, Moriarty JD. The effect of scaling titanium surfaces with metal and plastic instruments on cell attachment. *J Periodontol* 1990;61:491-496.
11. Fox SC, Moriarty JD, Kusy RP. The effects of scaling a titanium surface with metal and plastic instruments: An in vitro study. *J Periodontol* 1990;61:485-490.
12. Gantes BG, Nilveus R. The effects of different hygiene instruments on titanium surfaces: SEM observations. *Int J Periodont Rest Dent* 1991;11:225-239.
13. Homiak AW, Cook PA, DeBoer J. Effect of hygiene instrumentation on titanium abutments: A scanning electron microscopy study. *J Prosthet Dent* 1992;67:364-369.
14. Rapley JW, Swan RH, Hallmon WW, Millis MP. The surface characteristics produced by various oral hygiene instruments and materials on titanium implant abutments. *Int J Oral Maxillofac Implants* 1990;5:47-51.
15. Rühling A, Kocher T, Kreuzsch J, Plagman HC. Treatment of subgingival implant surfaces with Teflon coated sonic and ultrasonic scaler tips and various implant curettes. *Clin Oral Implants Res* 1994;5:19-29.
16. Thomson-Neal DM, Evans GH, Meffert RM. Effects of various prophylactic treatments on titanium sapphire and hydroxylapatite-coated implants: An SEM study. *Int J Periodont Rest Dent* 1989;9:301-311.
17. Parham PL, Cobb CM, French AA, Love JW, Drisko CL, Killoy WJ. Effects of an air-powder abrasive system on plasma-sprayed titanium implant surfaces: An in vitro evaluation. *J Oral Implantol* 1989;15:76-86.
18. Lie T, Leknes KN. Evaluation of the effect on root surfaces of air turbine scalers and ultrasonic instrumentation. *J Periodontol* 1985;56:522-531.